



**APEKS STERİLİZASYON**  
Medikal Hiz. San. ve Tic. Ltd. Şti.



# WHO WE ARE?

Apeks Sterilizasyon Medikal Hizmetleri Sanayi Ve Ticaret Ltd. Şti., which started off many years ago with the principle of being a pioneer in the sector, continues to work at full speed to offer you a more sterile life with its service quality and professional staff. Aiming to have an important share in the international market, our facility, which has an area of 800m<sup>2</sup>, has started to produce products with international standards, sterile and packaging services. Apeks family, which does not compromise professional service quality thanks to its harmony with developing technology, offers you a different experience.

As an innovator company, we are always seeking customer satisfaction. Encouraged by the customers, we added Surgical Clothing and Disposable Personal protective Equipment, surgical clothing to our production portfolio in 2022.



*Customer satisfaction*



*Professional staff*



*Delivery on time*



*Sterile packing*



*Roadway Connection  
Fast Sea Transport*

# PRODUCT LIST

---

- **DISPOSABLE SURGICAL DRAPES**

- **DISPOSABLE CLOTHINGS**

  - **PROTECTIVE COVERALLS**

  - Protective Coveralls Type 3/4
  - Protective Coveralls Type 5/6

## **SURGICAL GOWNS**

- Non Sterile Surgical Gowns
- Sterile Reinforced Surgical Gowns
- Sterile Standart Surgical Gowns
- Visitor Patient Gowns
- High Fluid Resistant Protective Gowns
- Scrub Suits

## **FFP2 MASKS**

- C Shape FFP2 Masks

## **SURGICAL MASKS**

- Earlooped Surgical Masks
- Tied Surgical Masks
- Pediatric Surgical Masks

## **THEATRE CAPS**

- Clip Caps
- Beret Style Theatre Caps
- Surgeon Hoods PE
- Surgeon Hoods PP
- Surgeon Hoods with Sweat Bands

## **CONSUMABLES**

- Laminated Boot Covers
- Laminated Non Slippery Boot Covers
- Shoe Covers PE
- Shoe Covers PP
- Overshoes PE
- Colonoscopy Shorts
- Gynecology Skirts
- Face Shields
- Bed Sheets
- Body Bags

# PROTECTIVE COVERALLS

APEKS1109-1110


Protective Coverall designed to help protect against Liquid Chemicals (Type 3-B, Type 4-B) hazardous dust (Type 5), certain light liquid splashes (Type 6)

Coveralls are commonly made of non-woven fabric that allows heat and sweat to leave the suit while preventing liquids and aerosols from entering it.

CE Marked

TYPE 3/4 COVERALLS PRODUCT CODE: 1110

TYPE 5/6 COVERALLS PRODUCT CODE: 1109




Biohazard



Liquid Protection



Liquid Protection



Dust Protection

Seamless shoulder and sleeve tops for fewer entry points for contaminants

Elasticized waist and ankles and knit or elastic cuffs on all models for increased comfort and mobility



# STERILIZATION MACHINE

Apeks Sterilizasyon Medikal Hizmetleri Sanayi Ve Ticaret Ltd. Şti., which started off many years ago with the principle of being a pioneer in the sector, continues to work at full speed to offer you a more sterile life with its service quality and professional staff. Aiming to have an important share in the international market, our facility, which has an area of 800m<sup>2</sup>, has started to produce products with international standards, sterile and packaging services. Apeks family, which does not compromise professional service quality thanks to its harmony with developing technology, offers you a different experience.



In our facility with 800m<sup>2</sup> working area, we serve as 100m<sup>2</sup> clean room, 400 m<sup>2</sup> EO Sterilization area, 100 m<sup>2</sup> ventilation area, 200m<sup>2</sup> packing and loading area



CE Marked  
ISO 13485:2016

## Clean Room

Packaging service is provided with sterilization rolls suitable for the size of the products. We have ISO 13485 certificate and it has been prepared in accordance with the conditions and rules of the document. Our daily product packaging capacity in the clean room is between 15,000 and 20,000

## E0 Sterilization area

We have 4 pcs 9m<sup>3</sup> and 1 pcs 14m<sup>3</sup> Ethylene Oxide Sterilization devices. Our sterilization team offers a daily capacity of 100,000-120,000 for 365 days. After the E0 Sterilization process, effective process control is provided by expert personnel with biological indicators.



## Ventilation area

It is the section where the products coming out of the sterilization device are ventilated for approximately 4 hours

## Packing and loading area

It is the area where the last process is carried out before the products reach the customer as a result of packaging, sterilization and ventilation. Your products are carefully packed and shipped to the customer



# SURGICAL DRAPE SETS

We provide surgical drape sets covering of medical procedures. Our drape sets include sterile and non-sterile surgical gowns, procedure suits, linens, covers, towels, and other accessories. We also offer individual surgical drape units and plain drapes (standard and customised).

General Surgery Drape Sets

Orthopaedic Surgery Drape Sets

Cardiology Surgery Drape Sets

Gynecology Surgery Drape Sets

Brain Surgery Drape Sets

Ophthalmology Drape Sets

Oncology Surgery Drape Sets

Urology Surgery Drape Sets

Radiology Surgery Drape Sets

E.N.T Surgery Drape Sets

Aesthetic Surgery Drape Sets

Dental Surgery Drape Sets



Anesthesia Drape



Chemotherapy Prep Mat



Drapes with Adhesive Tapes




# NON STERILE STANDARD SURGICAL GOWNS

Designed specifically for non-sterile procedures, this gown is constructed from SMS fabric and features a tie neck closure, knit cuffs, and a waist belt.

The materials used in the gowns is made of medical grade materials.

The knit cuffs of the gowns are thin and tight for ease of use. Breathable , non linting gowns provide comfortable usage to users.



CE Marked  
EN 13795



PRODUCT CODE : APEKS20 CODED PRODUCT SERIES





# STERILE STANDARD SURGICAL GOWNS

Sterile Standart Surgical Gowns are designed to use for short term operations, provides barrier between patient and surgeon. Sterilization prevents cross contamination.



Sterile Standart Surgical Gowns allows breathing and provides comfortable usage.

CE Marked

STERILE | EO



Various Material Thicknesses

PRODUCT CODE : APEKS21 CODED PRODUCT SERIES



# STERILE REINFORCED SURGICAL GOWNS

The materials used in the gowns is made of medical grade materials Gown fabric allows the skin to breathe, but it will remove alcohol, blood and similar body fluids from the environment and do not allow it to pass to the bottom.

They are deal for heavy procedures where strike-through and two-way contamination of the surgeon and patient could occur.

**CE Marked**    **STERILE | EO**



**PRODUCT CODE : APEKS22 CODED PRODUCT SERIES**

# PATIENT - VISITOR GOWNS

Patient-Visitor gowns provide solution against small amounts of liquid penetration and modesty to user

## Applicability

Basic care, Standard hospital medical unit



- One size fits all Single Use
- Tie Behind Neck Closure
- Latex-free

PRODUCT CODE : APEKS23 CODED PRODUCT SERIES

# LAMINATED PROTECTIVE GOWNS

Choosing the proper gowns for your surgical team members is an important part of your infection control strategy.

Impervious and High Fluid Gowns offer a high level of protection from blood and body fluids and help maintain the health and confidence of OR personnel.

Spunbond Fim Laminate with Polyethylene Coating provides more fluid resistance.

EN 14126



**Liquid  
Protection**



**Liquid  
Protection**



**Dust  
Protection**

PRODUCT CODE : APK20 CODED PRODUCT SERIES

# NON STERILE DISPOSABLE SCRUB SUIT

Always fresh and clean!

High cleanliness combined with a modern style for a fresh, hygienic start to every day. No wear and tear. No stains or odours. No need to worry.

CE Marked



40gsm Non Sterile medical grade material.

Elastic Cuffs

Without linting

Resistance to liquid penetration

With antistatic treatment

Breathable

Resistance to tensile strength

Ergonomic

Trousers strips closure system

All components are latex free.



PRODUCT CODE :2000 CODED PRODUCT SERIES

# RESPIRATORY PROTECTIVE MASKS

These products meet the requirements of EN 149: 2001 + A1: 2009 and EU Regulation 2016/425

Suitable Application Areas:

Hospitals/ Pharmaceutical Industry / Ceramics / Food / Asbestos / Hospital Laboratory / Virus Protection /Chemicals

Our masks are for single use, specially designed against collapse. Low breathing resistance filter technology ease breathing. Nose band helps to adjustment for each user.

Earlooped  
Latex-free  
CE Marked  
Regular Size



# SURGICAL MASKS

## TYPE II R

Comfortable, Latex-free

Adaptable nose bar for better fitting.

Microbially cleaned by UV before packing. Masks does not contain allergic substance. Provides easy breathing.

Thermosealed



REF:1717A



REF:1717



REF:1717K

# SURGICAL CAPS

Surgical Cap helps protect the healthcare facility and healthcare worker from microorganism contamination.




**ASTRO CAPS-APEKS1002**



**SURGEON HOODS PE-APEKS1003**



**CLIP CAPS-APEKS1000**



**BERET STYLE CAPS -APEKS1001**



# CONSUMABLES

---



**Laminated Booth Covers**



**Shoe Covers**



**Face Shield**



**Gynecology Skirt**



**Colonoscopy Short**



**Bed Sheet**



**Body Bag**



**APEKS STERİLİZASYON**  
Medikal Hiz. San. ve Tic. Ltd. Şti.

**APEKS STERİLİZASYON MEDİKAL HİZMETLERİ SANAYİ VE TİCARET LİMİTED ŞİRKETİ**

**Tel: 0(216) 540 4009**

**Şerifali Mahallesi, Serdivan Sokak, No:28 Ümraniye/İSTANBUL**

**[www.apekssterilizasyon.com](http://www.apekssterilizasyon.com)**

**[info@apekssterilizasyon.com](mailto:info@apekssterilizasyon.com)**



**ISO 13485:2016**



**ISO 9001:2015**



**ISO 45001:2018**